


FILM LAYER: 5.6 MIL THICK, HIGH CLARITY, HYDROPHOBIC

ADHESIVE LAYER: 1 MIL THICK, PRESSURE SENSITIVE ACRYLIC

RELEASE LAYER: 90 LB. POLY COATED KRAFT

PACKAGED: 50/BAG



**Art Robbins
Instruments**

TITLE
SEAL
HANGING DROP

SIZE	CAGE CODE	DWG NO	REV
B		600-4006-00	A
SCALE	SHEET		
A			