



# INTELLI-PLATE®


## 96-well INTELLI-PLATE®




Catalog #: 102-0011-00  
Original INTELLI-PLATE  
Well Volume: 4µL and 10µL  
Reservoir Volume: 300µL max




Catalog #: 102-0001-00 \*\*  
96-2 Low volume reservoir  
Well Volume: 4µL and 10µL  
Reservoir Volume: 100µL max




Catalog #: 102-0001-10\*\*  
96-2 Low profile  
Well Volume: 4µL and 10µL  
Reservoir Volume: 100µL



Catalog #: 102-0001-03\*\*  
96-3 Low volume reservoir  
Well Volume: 1µL  
Reservoir Volume: 100µL max




Catalog #: 102-0001-13\*\*  
96-3 low profile INTELLI-PLATE  
Well Volume: 1µL  
Reservoir Volume: 100µL max




Catalog #: 102-0001-01  
Flat bottom  
Well Volume: N/A  
Reservoir Volume: 300µL max

\*\* UV Compatible Plate




Catalog #: 102-0001-20\*\*  
96-2 Shallow well  
Well Volume: 2µL  
Reservoir Volume: 100µL max




Catalog #: 102-0001-21\*\*  
96-2 Shallow well low profile  
Well Volume: 2µL  
Reservoir Volume: 100µL max


## 24-well INTELLI-PLATE® 48-well INTELLI-PLATE®



Catalog #: 102-0004-00  
24-4  
Well Volume: 5µL  
Reservoir Volume: 650µL max



Catalog #: 102-0003-00  
48-3 Screening  
Well Volume: 4µL  
Reservoir Volume: 500µL max



Catalog #: 102-0002-00  
48-2 Optimization  
Well Volume: 4µL and 20µL  
Reservoir Volume: 500µL max

INTELLI-PLATE is a registered trademark of Art Robbins Instruments, LLC.



For more information, contact us at: Tel: (888)658-5300

[www.artrobbins.com](http://www.artrobbins.com)

[info@artrobbins.com](mailto:info@artrobbins.com)